

For Research Use Only

KLHL4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84392-2-PBS



Basic Information

Catalog Number:

84392-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15529

GenBank Accession Number:

BC142652

GeneID (NCBI):

56062

UNIPROT ID:

Q9C0H6

Full Name:

kelch-like 4 (Drosophila)

Calculated MW:

720 aa, 81 kDa

Purification Method:

Protein A purification

CloneNo.:

241773E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01263-1, KLHL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84392-2-PBS. Detection antibody: 84392-1-PBS. Standard: Ag15529. Range: 0.781-100 ng/mL.