

For Research Use Only

# Mouse IFNGR1/CD119 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84415-6-PBS



## Basic Information

Catalog Number:	84415-6-PBS	GenBank Accession Number:	NM_010511.3	Purification Method:	Protein A purification
Concentration:	1 mg/ml	GeneID (NCBI):	15979	CloneNo.:	241832H4
Source:	Rabbit	UNIPROT ID:	P15261		
Isotype:	IgG	Full Name:	interferon gamma receptor 1		
Immunogen Catalog Number:	EG1536	Calculated MW:	52kDa		

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
mouse

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01291-3, MOUSE IFN-gamma R1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84415-2-PBS. Detection antibody: 84415-6-PBS. Standard: Eg1536. Range: 0.625-80 ng/mL