

For Research Use Only

# Rat NKp46/NCR1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84416-1-PBS



## Basic Information

Catalog Number:

84416-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_057199.1

GeneID (NCBI):

117547

UNIPROT ID:

Q9Z0H5

Full Name:

natural cytotoxicity triggering  
receptor 1

Calculated MW:

37kDa

Purification Method:

Protein A purification

CloneNo.:

241811C4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01299-3, RAT NKp46/NCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84416-3-PBS. Detection antibody: 84416-1-PBS. Standard: Eg1469. Range: 0.313-40 ng/mL