

For Research Use Only

Rat IL-4 Recombinant antibody, PBS Only

Catalog Number: 84422-4-PBS



Basic Information

Catalog Number:

84422-4-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1113

GenBank Accession Number:

NM_201270.1

GeneID (NCBI):

287287

UNIPROT ID:

P20096

Full Name:

interleukin 4

Calculated MW:

16kDa

Purification Method:

Protein A purification

CloneNo.:

241249G7

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

