

For Research Use Only

Mouse TWEAK R/TNFRSF12 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84423-1-PBS



Basic Information

Catalog Number:

84423-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_013749.2

GeneID (NCBI):

27279

UNIPROT ID:

Q9CR75

Full Name:

tumor necrosis factor receptor
superfamily, member 12a

Calculated MW:

14kDa

Purification Method:

Protein A purification

CloneNo.:

241646C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01283-3, MOUSE TWEAK R/TNFRSF12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84423-1-PBS. Detection antibody: 84423-3-PBS. Standard: Eg1838. Range: 0.156-20 ng/mL