

For Research Use Only

# PSMA5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84424-2-PBS



## Basic Information

Catalog Number:

84424-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34928

GenBank Accession Number:

BC103751

GeneID (NCBI):

5686

UNIPROT ID:

P28066

Full Name:

proteasome (prosome, macropain)

subunit, alpha type, 5

Calculated MW:

241 aa, 26 kDa

Purification Method:

Protein A purification

CloneNo.:

241629E8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01286-1, PSMA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84424-2-PBS. Detection antibody: 84424-1-PBS. Standard: Ag34928. Range: 0.156-20 ng/mL.



Cytometric bead array standard curve of MP01286-2, PSMA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84424-2-PBS. Detection antibody: 84424-3-PBS. Standard: Ag34928. Range: 0.156-20 ng/mL.