

For Research Use Only

HSPA13 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 84426-3-PBS



Basic Information

Catalog Number:

84426-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3349

GenBank Accession Number:

BC036370

GeneID (NCBI):

6782

UNIPROT ID:

P48723

Full Name:

heat shock protein 70kDa family,
member 13

Calculated MW:

471 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

241572E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean HSPA13 concentration was determined to be 24.8 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 13.9 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 6.6 ng/mL in K-562 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP01302-3, HSPA13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84426-3-PBS. Detection antibody: 84426-1-PBS. Standard: Ag3349. Range: 0.781-100 ng/mL.