

For Research Use Only

CXCL9/MIG Recombinant antibody, PBS Only (Detector)

Catalog Number: 84427-1-PBS



Basic Information

Catalog Number:	84427-1-PBS	GenBank Accession Number:	BC095396	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	4283	CloneNo.:	241437A2
Source:	Rabbit	UNIPROT ID:	Q07325		
Isotype:	IgG	Full Name:	chemokine (C-X-C motif) ligand 9		
Immunogen Catalog Number:	AG17932	Calculated MW:	14 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01308-1, CXCL9/MIG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84427-3-PBS. Detection antibody: 84427-1-PBS. Standard: Ag17932. Range: 0.781-100 ng/mL.