

For Research Use Only

# Mouse SIRP alpha/CD172a Recombinant antibody, PBS Only (Capture)

Catalog Number: 84430-4-PBS



## Basic Information

Catalog Number:

84430-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_007547.4

GeneID (NCBI):

19261

UNIPROT ID:

P97797-2

Full Name:

signal-regulatory protein alpha

Calculated MW:

56 kDa

Purification Method:

Protein A purification

CloneNo.:

241794G8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01311-1, MOUSE SIRP Alpha/CD172a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84430-4-PBS. Detection antibody: 84430-3-PBS. Standard: Eg1520. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01311-3, MOUSE SIRP Alpha/CD172a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84430-4-PBS. Detection antibody: 84430-1-PBS. Standard: Eg1520. Range: 0.781-100 ng/mL