

For Research Use Only

# Mouse Ki-67 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84432-2-PBS



## Basic Information

Catalog Number:

84432-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001081117.2

GeneID (NCBI):

17345

UNIPROT ID:

E9PVX6-1

Full Name:

antigen identified by monoclonal  
antibody Ki 67

Calculated MW:

351kDa

Purification Method:

Protein A purification

CloneNo.:

241768H7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01277-1, MOUSE Ki-67 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84432-2-PBS. Detection antibody: 84432-1-PBS. Standard: Eg1548. Range: 0.625-40 ng/mL.