

For Research Use Only

Mouse 5'-Nucleotidase/CD73 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84434-2-PBS



Basic Information

Catalog Number:

84434-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_011851.4

GeneID (NCBI):

23959

UNIPROT ID:

Q61503

Full Name:

5' nucleotidase, ecto

Calculated MW:

64kDa

Purification Method:

Protein A purification

CloneNo.:

241541A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01279-2, MOUSE CD73 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84434-5-PBS. Detection antibody: 84434-2-PBS. Standard: Eg1389. Range: 0.625-40 ng/mL.