

For Research Use Only

# DDR1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84438-2-PBS



## Basic Information

Catalog Number:

84438-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001297654.2

GeneID (NCBI):

780

UNIPROT ID:

Q08345-1

Full Name:

discoidin domain receptor tyrosine kinase 1

Calculated MW:

101kDa

Purification Method:

Protein A purification

CloneNo.:

241796C6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01290-1, DDR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84438-2-PBS. Detection antibody: 84438-3-PBS. Standard:null. Range: null.



Cytometric bead array standard curve of MP01290-2, DDR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84438-2-PBS. Detection antibody: 84438-1-PBS. Standard: Eg1675. Range: 0.313-40 ng/mL.