

For Research Use Only

Mouse Nectin-2/CD112 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84439-1-PBS



Basic Information

Catalog Number:

84439-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2019

GenBank Accession Number:

NM_008990.3

GeneID (NCBI):

19294

UNIPROT ID:

P32507-1

Full Name:

poliovirus receptor-related 2

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

241798A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Nectin-2/CD112 concentration was determined to be 9.6 ng/mL in RAW264.7 cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP01292-1, MOUSE Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84439-3-PBS. Detection antibody: 84439-1-PBS. Standard: Eg2019. Range: 0.313-40 ng/mL.