

For Research Use Only

Nectin-2/CD112 Recombinant antibody

Catalog Number: 84439-1-RR



Basic Information

Catalog Number: 84439-1-RR	GenBank Accession Number: NM_008990.3	Purification Method: Protein A purification
Size: Source: Rabbit	GeneID (NCBI): 19294	CloneNo.: 241798A2
Isotype: IgG	UNIPROT ID: P32507-1	
	Full Name: poliovirus receptor-related 2	
	Calculated MW: 57kDa	

Applications

Tested Applications:
ELISA
Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data