

For Research Use Only

CDH15 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84440-2-PBS



Basic Information

Catalog Number:

84440-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25763

GenBank Accession Number:

BC008951

GeneID (NCBI):

1013

UNIPROT ID:

P55291

Full Name:

cadherin 15, type 1, M-cadherin (myotubule)

Calculated MW:

89 kDa

Purification Method:

Protein A purification

CloneNo.:

241837B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01295-2, M-cadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84440-2-PBS. Detection antibody: 84440-1-PBS. Standard: Ag25763. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01295-3, M-cadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84440-4-PBS. Detection antibody: 84440-2-PBS. Standard: Ag25763. Range: 0.625-40 ng/mL