

For Research Use Only

# P2RX6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84458-2-PBS



## Basic Information

Catalog Number:

84458-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12331

GenBank Accession Number:

BC033488

GeneID (NCBI):

9127

UNIPROT ID:

O15547

Full Name:

purinergic receptor P2X, ligand-gated ion channel, 6

Calculated MW:

431 aa, 48 kDa

Purification Method:

Protein A purification

CloneNo.:

241620C6

## Applications

Tested Applications:

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01327-1, P2RX6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84458-2-PBS. Detection antibody: 84458-1-PBS. Standard: Ag12331. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01327-2, P2RX6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84458-2-PBS. Detection antibody: 84458-3-PBS. Standard: Ag12331. Range: 0.781-100 ng/mL.