

For Research Use Only

DC-SIGNR/CD299 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84461-3-PBS



Basic Information

Catalog Number:

84461-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC038851

GeneID (NCBI):

10332

UNIPROT ID:

Q9H2X3

Full Name:

C-type lectin domain family 4, member M

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

241682H5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01330-1, DC-SIGNR/CD299 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84461-2-PBS. Detection antibody: 84461-3-PBS. Standard: RP00011. Range: 0.313-40 ng/mL.