

For Research Use Only

CD100 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84465-2-PBS



Basic Information

Catalog Number:

84465-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC054500

GeneID (NCBI):

10507

UNIPROT ID:

Q92854

Full Name:

sema domain, immunoglobulin domain (Ig), transmembrane domain (TM) and short cytoplasmic domain, (semaphorin) 4D

Calculated MW:

96 kDa

Purification Method:

Protein A purification

CloneNo.:

241709C5

Applications

Tested Applications:

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01337-2, CD100 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84465-2-PBS. Detection antibody: 84465-3-PBS. Standard: RP02354. Range: 0.625-80 ng/mL