

For Research Use Only

ADAMTS2 Recombinant antibody, PBS Only

Catalog Number: 84476-6-PBS



Basic Information

Catalog Number:

84476-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34957

GenBank Accession Number:

NM_014244

GeneID (NCBI):

9509

UNIPROT ID:

O95450

Full Name:

ADAM metalloproteinase with
thrombospondin type 1 motif, 2

Calculated MW:

135 kDa

Observed MW:

110-130 kDa

Purification Method:

Protein A purification

CloneNo.:

241860B8

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

A disintegrin and metalloproteinase with thrombospondin motifs 2 (ADAMTS-2) is a metalloproteinase that plays a key role in the processing of fibrillar procollagen precursors into mature collagen molecules by excising the amino-propeptide. ADAMTS-2 can, like the procollagen C-proteinases, be regulated by transforming growth factor- β 1 (TGF- β 1), with implications for mechanisms whereby this growth factor effects net increases in formation of extracellular matrix. Mutations in the ADAMTS-2 gene cause dermatosparaxis in human (dermatosparactic type of Ehlers-Danlos syndrome, also known as Ehlers-Danlos syndrome type VIIC). Human ADAMTS-2 has eight potential sites for Asn-linked glycosylation and other mammals (PMID: 16046392, 12646579).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84476-6-RR (ADAMTS2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84476-6-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A431 cells using ADAMTS2 antibody (84476-6-RR, Clone: 241860B8) at dilution of 1:250 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 84476-6-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84476-6-RR against Human ADAMTS2 were performed. The affinity constant is 0.472 nM.