

For Research Use Only

CCDC25 Recombinant antibody

Catalog Number: 84484-3-RR



Basic Information

Catalog Number:

84484-3-RR

Size:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15593

GenBank Accession Number:

BC032588

GeneID (NCBI):

55246

UNIPROT ID:

Q86WRO

Full Name:

coiled-coil domain containing 25

Calculated MW:

208 aa, 24 kDa

Observed MW:

24-26 kDa

Purification Method:

Protein A purification

CloneNo.:

241878H5

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : A549 cells, HeLa cells, Jurkat cells, A2780 cells, K-562 cells, MCF-7 cells, RAW 264.7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84484-3-RR (CCDC25 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 84484-3-RR against Human CCDC25 were performed. The affinity constant is 0.771 nM.