

For Research Use Only

# GINS2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84486-1-PBS



## Basic Information

Catalog Number:

84486-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9180

GenBank Accession Number:

BC010164

GeneID (NCBI):

51659

UNIPROT ID:

Q9Y248

Full Name:

GINS complex subunit 2 (Psf2 homolog)

Calculated MW:

185 aa, 21 kDa

Purification Method:

Protein A purification

CloneNo.:

241870B11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01325-1, GINS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84486-2-PBS. Detection antibody: 84486-1-PBS. Standard: Ag9180. Range: 0.313-20 ng/mL