

For Research Use Only

# Mouse EMMPRIN/CD147 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84488-2-PBS



## Basic Information

Catalog Number:

84488-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_009768.2

GeneID (NCBI):

12215

UNIPROT ID:

P18572-1

Full Name:

basigin

Calculated MW:

42kDa

Purification Method:

Protein A purification

CloneNo.:

241852B9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01320-1, MOUSE CD147 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84488-2-PBS. Detection antibody: 84488-4-PBS. Standard: Eg1530. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP01320-2, MOUSE CD147 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84488-1-PBS. Detection antibody: 84488-2-PBS. Standard: Eg1530. Range: 0.781-50 ng/mL