

For Research Use Only

UHRF1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84489-3-PBS



Basic Information

Catalog Number:

84489-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17498

GenBank Accession Number:

BC113875

GeneID (NCBI):

29128

UNIPROT ID:

Q96T88

Full Name:

ubiquitin-like with PHD and ring
finger domains 1

Calculated MW:

806 aa, 91 kDa

Purification Method:

Protein A purification

CloneNo.:

241872G2

Applications

Tested Applications:

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01341-2, UHRF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84489-1-PBS. Detection antibody: 84489-3-PBS. Standard: Ag17498. Range: 0.156-20 ng/mL