

For Research Use Only

Mouse CXCR4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84506-2-PBS



Basic Information

Catalog Number:

84506-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009911.3

GeneID (NCBI):

12767

UNIPROT ID:

P70658

Full Name:

chemokine (C-X-C motif) receptor 4

Calculated MW:

40kDa

Purification Method:

Protein A purification

CloneNo.:

241956H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01396-2, MOUSE CXCR4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84506-1-PBS. Detection antibody: 84506-2-PBS. Standard: Eg2231. Range: 0.156-20 ng/mL.