

For Research Use Only

CXCR4 Recombinant antibody

Catalog Number: 84506-2-RR



Basic Information

Catalog Number:

84506-2-RR

Size:

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009911.3

GeneID (NCBI):

12767

UNIPROT ID:

P70658

Full Name:

chemokine (C-X-C motif) receptor 4

Calculated MW:

40kDa

Purification Method:

Protein A purification

CloneNo.:

241956H2

Applications

Tested Applications:

ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

