

For Research Use Only

CD37 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84509-3-PBS



Basic Information

Catalog Number:

84509-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1764

GenBank Accession Number:

NM_001774.3

GeneID (NCBI):

951

UNIPROT ID:

P11049-1

Full Name:

CD37 molecule

Calculated MW:

32kDa

Purification Method:

Protein A purification

CloneNo.:

242011B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01362-1, CD37 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84509-3-PBS. Detection antibody: 84509-2-PBS. Standard: Eg1764. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01362-2, CD37 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84509-3-PBS. Detection antibody: 84509-1-PBS. Standard: Eg1764. Range: 0.313-40 ng/mL