

For Research Use Only

# Mouse IL-3RB/CD131 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84513-2-PBS



## Basic Information

|  |  |   |
|--|--|---|
| <b>Catalog Number:</b><br>84513-2-PBS      | <b>GenBank Accession Number:</b><br>NM_007780.4  | <b>Purification Method:</b><br>Protein A purification |
| <b>Concentration:</b><br>1 mg/ml           | <b>GeneID (NCBI):</b><br>12983   | <b>CloneNo.:</b><br>241954C11                         |
| <b>Source:</b><br>Rabbit                   | <b>UNIPROT ID:</b><br>P26955   |   |
| <b>Isotype:</b><br>IgG                     | <b>Full Name:</b><br>colony stimulating factor 2 receptor, beta, low-affinity (granulocyte-macrophage) |   |
| <b>Immunogen Catalog Number:</b><br>EG2173 | <b>Calculated MW:</b><br>99kDa   |   |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
mouse

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01393-2, MOUSE IL-3RB/CD131 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84513-2-PBS. Detection antibody: 84513-3-PBS. Standard: Eg2173. Range: 0.313-40 ng/mL