

For Research Use Only

Mouse IL-3RB/CD131 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84513-3-PBS



Basic Information

Catalog Number: 84513-3-PBS	GenBank Accession Number: NM_007780.4	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 12983	CloneNo.: 241954811
Source: Rabbit	UNIPROT ID: P26955	
Isotype: IgG	Full Name: colony stimulating factor 2 receptor, beta, low-affinity (granulocyte-macrophage)	
Immunogen Catalog Number: EG2173	Calculated MW: 99kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01393-1, MOUSE IL-3RB/CD131 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84513-1-PBS. Detection antibody: 84513-3-PBS. Standard: Eg2173. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01393-2, MOUSE IL-3RB/CD131 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84513-2-PBS. Detection antibody: 84513-3-PBS. Standard: Eg2173. Range: 0.313-40 ng/mL