

For Research Use Only

# DCIR/CLEC4A Recombinant antibody, PBS Only (Capture)

Catalog Number: 84515-4-PBS



## Basic Information

|                                       |   |   |
|---------------------------------------|---|---|
| <b>Catalog Number:</b><br>84515-4-PBS | <b>GenBank Accession Number:</b><br>NM_016184.4                 | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1 mg/ml               | <b>GeneID (NCBI):</b><br>50856                                  | <b>CloneNo.:</b><br>241854G1                          |
| <b>Source:</b><br>Rabbit              | <b>UNIPROT ID:</b><br>Q9UMR7-1                                  |   |
| <b>Isotype:</b><br>IgG                | <b>Full Name:</b><br>C-type lectin domain family 4,<br>member A |   |
|                                       | <b>Calculated MW:</b><br>28kDa                                  |   |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01361-1, CLEC4A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84515-4-PBS. Detection antibody: 84515-1-PBS. Standard: Eg1350. Range: 0.469-60 ng/mL.