

For Research Use Only

TRPC6 Recombinant monoclonal antibody, PBS Only

Catalog Number: 84517-5-PBS



Basic Information

Catalog Number:

84517-5-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC093658

GeneID (NCBI):

7225

UNIPROT ID:

Q9Y210

Full Name:

transient receptor potential cation channel, subfamily C, member 6

Calculated MW:

931 aa, 106 kDa

Observed MW:

130 kDa

Purification Method:

Protein A purification

CloneNo.:

251552A8

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Transient receptor potential canonical (TRPC) channels mediate the influx of different types of cations through the cell membrane and are involved in many functions of the organism (PMID: 28508319). The TRPC subfamily consists of seven members, TRPC1, TRPC2, TRPC3, TRPC4, TRPC5, TRPC6, TRPC7 (PMID: 32872338). TRPC6 is a Ca(2+)-permeable non-selective cation channel expressed in most human tissues. TRPC6 has been considered as potential therapeutic targets for malignant glioma treatment (PMID: 35489424).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84517-5-RR (TRPC6 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84517-5-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A549 cells using TRPC6 antibody (84517-5-RR, Clone: 251552A8) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 84517-5-PBS in a different storage buffer formulation.