

For Research Use Only

Mouse VE-cadherin Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84529-2-PBS



Basic Information

Catalog Number:

84529-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009868.4

GeneID (NCBI):

12562

UNIPROT ID:

P55284

Full Name:

cadherin 5

Calculated MW:

88kDa

Purification Method:

Protein A purification

CloneNo.:

241953C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01386-1, MOUSE VE-cadherin/CD144 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84529-2-PBS. Detection antibody: 84529-1-PBS. Standard: Eg0921. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01386-2, MOUSE VE-cadherin/CD144 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84529-3-PBS. Detection antibody: 84529-2-PBS. Standard: Eg0921. Range: 0.781-100 ng/mL