

For Research Use Only

IL-1RA Recombinant antibody, PBS Only (Capture)

Catalog Number: 84532-1-PBS



Basic Information

Catalog Number:	84532-1-PBS	GenBank Accession Number:	BC009745	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	3557	CloneNo.:	241879D3
Source:	Rabbit	UNIPROT ID:	P18510		
Isotype:	IgG	Full Name:	interleukin 1 receptor antagonist		
Immunogen Catalog Number:	AG1277	Calculated MW:	20 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01389-2, IL-1RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84532-1-PBS. Detection antibody: 84532-3-PBS. Standard: Ag1277. Range: 0.625-80 ng/mL