

For Research Use Only

IL17RA Recombinant antibody, PBS Only (Capture)

Catalog Number: 84535-3-PBS



Basic Information

Catalog Number:

84535-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC011624

GeneID (NCBI):

23765

UNIPROT ID:

Q96F46

Full Name:

interleukin 17 receptor A

Calculated MW:

866 aa, 96 kDa

Purification Method:

Protein A purification

CloneNo.:

241891G7

Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA,
Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP01397-3, Human IL17RA Recombinant Matched Antibody Pair - PBS only. 84535-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2317. 84535-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



The mean IL17RA concentration was determined to be 8.86 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP01397-2, IL17RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84535-3-PBS. Detection antibody: 84535-1-PBS. Standard: Eg2317. Range: 0.625-80 ng/mL