

For Research Use Only

FGF19 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84536-4-PBS



Basic Information

Catalog Number:

84536-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

HZ-1330

GenBank Accession Number:

BC017664

GeneID (NCBI):

9965

UNIPROT ID:

O95750

Full Name:

fibroblast growth factor 19

Calculated MW:

216 aa, 24 kDa

Purification Method:

Protein A purification

CloneNo.:

241740E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01398-2, FGF19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84536-1-PBS. Detection antibody: 84536-4-PBS. Standard: HZ-1330. Range: 0.078-20 ng/mL