

For Research Use Only

# Mouse F4/80 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84538-3-PBS



## Basic Information

Catalog Number:

84538-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_010130.4

GeneID (NCBI):

13733

UNIPROT ID:

Q61549

Full Name:

EGF-like module containing, mucin-like, hormone receptor-like sequence  
1

Calculated MW:

102 kDa

Purification Method:

Protein A purification

CloneNo.:

241959F9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01359-1, MOUSE F4/80 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84538-3-PBS. Detection antibody: 84538-2-PBS. Standard: Eg0788. Range: 0.625-40 ng/mL.