

For Research Use Only

BAFF Recombinant antibody

Catalog Number: 84541-3-RR



Basic Information

Catalog Number:

84541-3-RR

Size:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_033622

GeneID (NCBI):

24099

UNIPROT ID:

Q9WU72

Full Name:

tumor necrosis factor (ligand)
superfamily, member 13b

Calculated MW:

34kd

Observed MW:

32kd, 17kd

Purification Method:

Protein A purification

CloneNo.:

242006E4

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : U-937 cells, Jurkat cells, RAW 264.7 cells

Background Information

BAFF (B-cell-activating factor), also called Tnfsf13b (Tumor necrosis factor ligand superfamily member 13B) and CD257, is a cytokine that binds to TNFRSF13B/TACI and TNFRSF17/BCMA. TNFSF13/APRIL binds to the same 2 receptors. Together, they form a 2 ligands -2 receptors pathway involved in stimulating B- and T-cell function and the regulation of humoral immunity. A third B-cell-specific BAFF-receptor (BAFFR/BR3) promotes the survival of mature B-cells and the B-cell response (PMID: 12867412). BAFF cleaved into 2 chains: membrane form is 32-35 kDa; soluble form is 17 kDa.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84541-3-RR (BAFF antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



U-937 cells were subjected to SDS PAGE followed by western blot with 84541-3-RR (BAFF antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.