

For Research Use Only

ER Recombinant antibody, PBS Only

Catalog Number: 84564-1-PBS



Basic Information

Catalog Number:

84564-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2234

GenBank Accession Number:

BC128573

GeneID (NCBI):

2099

UNIPROT ID:

P03372

Full Name:

estrogen receptor 1

Calculated MW:

595 aa, 66 kDa

Purification Method:

Protein A purification

CloneNo.:

241825C1

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assay of 84564-1-PBS against Human ER was performed. The affinity constant is 2.62 nM.