

For Research Use Only

DPCR1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84567-1-PBS



Basic Information

Catalog Number:

84567-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33318

GenBank Accession Number:

NM_080870

GeneID (NCBI):

135656

UNIPROT ID:

Q3MIW9

Full Name:

diffuse panbronchiolitis critical
region 1

Calculated MW:

56 kDa

Purification Method:

Protein A purification

CloneNo.:

241974B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01409-1, DPCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84567-1-PBS. Detection antibody: 84567-3-PBS. Standard: Ag33318. Range: 1.56-100 ng/mL