

For Research Use Only

TRAF2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84581-2-PBS



Basic Information

Catalog Number:

84581-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25387

GenBank Accession Number:

BC032410

GeneID (NCBI):

7186

UNIPROT ID:

Q12933

Full Name:

TNF receptor-associated factor 2

Calculated MW:

56 kDa

Purification Method:

Protein A purification

CloneNo.:

241899G11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01404-1, TRAF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84581-2-PBS. Detection antibody: 84581-1-PBS. Standard: Ag25387. Range: 1.563-200 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84581-2-PBS against Human TRAF2 was performed. The affinity constant is below 1 pM.