

For Research Use Only

IL15RA Recombinant antibody, PBS Only (Detector)

Catalog Number: 84582-2-PBS



Basic Information

Catalog Number:

84582-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC107777

GeneID (NCBI):

3601

UNIPROT ID:

Q13261

Full Name:

interleukin 15 receptor, alpha

Calculated MW:

231 aa, 24 kDa

Purification Method:

Protein A purification

CloneNo.:

241756A7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01405-1, IL-15RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84582-4-PBS. Detection antibody: 84582-2-PBS. Standard: RP00015. Range: 0.469-60 ng/mL