

For Research Use Only

# Granzyme A Recombinant antibody, PBS Only (Detector)

Catalog Number: 84583-1-PBS



## Basic Information

Catalog Number:

84583-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1818

GenBank Accession Number:

BC015739

GeneID (NCBI):

3001

UNIPROT ID:

P12544

Full Name:

granzyme A (granzyme 1, cytotoxic T-lymphocyte-associated serine esterase 3)

Calculated MW:

29 kDa

Purification Method:

Protein A purification

CloneNo.:

241979B5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01417-1, Granzyme A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84583-2-PBS. Detection antibody: 84583-1-PBS. Standard: Ag1818. Range: 0.625-80 ng/mL