

For Research Use Only

BCAP31 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 84584-4-PBS



Basic Information

Catalog Number:

84584-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1668

GenBank Accession Number:

BC014323

GeneID (NCBI):

10134

UNIPROT ID:

P51572

Full Name:

B-cell receptor-associated protein 31

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

241978H5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01418-2, BAP31 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84584-3-PBS. Detection antibody: 84584-4-PBS. Standard: Ag1668. Range: 0.625-80 ng/mL



Biolayer interferometry (BLI) kinetic assay of 84584-4-PBS against Human BCAP31 was performed. The affinity constant is below 1 pM.