

For Research Use Only

CXCR6 Recombinant antibody, PBS Only

Catalog Number: 84591-1-PBS



Basic Information

Catalog Number:

84591-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC033584

GeneID (NCBI):

10663

UNIPROT ID:

O00574

Full Name:

chemokine (C-X-C motif) receptor 6

Calculated MW:

39 kDa

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

241662E5

Applications

Tested Applications:

WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

CXCR6, or C-X-C chemokine receptor 6, is a G protein-coupled receptor with seven transmembrane domains that belong to the CXCR chemokine receptor family. This family also includes CXCR1, CXCR2, CXCR3, CXCR4, CXCR5, and CXCR7. CXCR6 is expressed in several T lymphocyte subsets and bone marrow stromal cells. The encoded protein and its exclusive ligand, chemokine ligand 16 (CCL16), are part of a signaling pathway that regulates T lymphocyte migration to various peripheral tissues such as the liver, spleen red pulp, intestine, lungs, and skin. CXCR6 plays a crucial role in immunosurveillance through its interaction with epithelial cells and certain subpopulations of intratumor dendritic cells (DC), particularly DC3, which express CXCL16 at the membrane.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84591-1-RR (CXCR6 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84591-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 84591-1-RR (CXCR6 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 84591-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of un-fixed human lung cancer tissue using CXCR6 antibody (84591-1-RR, Clone: 241662E5) at dilution of 1:200 and Multi-rAb CoraLite ® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002). This data was developed using the same antibody clone with 84591-1-PBS in a different storage buffer formulation.