

For Research Use Only

CD68 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 84596-9-PBS



Basic Information

Catalog Number:

84596-9-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2061

GenBank Accession Number:

BC015557

GeneID (NCBI):

968

UNIPROT ID:

P34810

Full Name:

CD68 molecule

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

261327C10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD68 concentration was determined to be 6.1 ng/mL in human serum.



Cytometric bead array standard curve of MP01429-4, CD68 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84596-10-PBS. Detection antibody: 84596-9-PBS. Standard: Eg2061. Range: 0.625-80 ng/mL.