

For Research Use Only

Mouse PDGF R alpha/PDGFR Recombinant antibody, PBS Only (Capture)

Catalog Number: 84603-2-PBS



Basic Information

Catalog Number:

84603-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001083316.2

GeneID (NCBI):

18595

UNIPROT ID:

P26618-1

Full Name:

platelet derived growth factor
receptor, alpha polypeptide

Calculated MW:

123kDa

Purification Method:

Protein A purification

CloneNo.:

242025E4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01458-1, MOUSE PDGF R alpha/PDGFRA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84603-2-PBS. Detection antibody: 84603-1-PBS. Standard: Eg1625. Range: 0.469-60 ng/mL