

For Research Use Only

Pregnancy Zone Protein Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84605-1-PBS



Basic Information

Catalog Number:

84605-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16471

GenBank Accession Number:

BC111756

GeneID (NCBI):

5858

Full Name:

pregnancy-zone protein

Calculated MW:

1482 aa, 164 kDa

Purification Method:

Protein A purification

CloneNo.:

242114C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01425-1, Pregnancy Zone Protein Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84605-1-PBS. Detection antibody: 84605-3-PBS. Standard: Ag16471. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01425-2, Pregnancy Zone Protein Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84605-2-PBS. Detection antibody: 84605-1-PBS. Standard: Ag16471. Range: 0.781-100 ng/mL