

For Research Use Only

EMR2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84608-4-PBS



Basic Information

Catalog Number:

84608-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_013447.3

GeneID (NCBI):

30817

UNIPROT ID:

Q9UHX3-1

Full Name:

egf-like module containing, mucin-like, hormone receptor-like 2

Calculated MW:

90 kDa

Purification Method:

Protein A purification

CloneNo.:

242042G9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01442-3, Human EMR2 Recombinant Matched Antibody Pair - PBS only. 84608-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2479. 84608-4-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.



The mean EMR2 concentration was determined to be 93.8 ng/mL in HL-60 cell extract based on a 1.4 mg/mL extract load, 98.6 ng/mL in THP-1 cell extract based on a 4.4 mg/mL extract load, 27.8 ng/mL in K-562 cell extract based on a 2.4 mg/mL extract load and 28.5 ng/mL in Jurkat cell extract based on a 4.2 mg/mL extract load.