

For Research Use Only

TTC9C Recombinant antibody, PBS Only (Detector)

Catalog Number: 84611-1-PBS



Basic Information

Catalog Number:

84611-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37081

GenBank Accession Number:

BC032123

GeneID (NCBI):

283237

UNIPROT ID:

Q8N5M4

Full Name:

tetratricopeptide repeat domain 9C

Purification Method:

Protein A purification

CloneNo.:

242061C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01448-1, TTC9C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84611-2-PBS. Detection antibody: 84611-1-PBS. Standard: Ag37081. Range: 0.781-100 ng/mL.