

For Research Use Only

Mouse Fibrinogen gamma chain Recombinant antibody, PBS Only (Detector)

Catalog Number: 84616-3-PBS



Basic Information

Catalog Number:

84616-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1837

GenBank Accession Number:

NM_133862.2

GeneID (NCBI):

99571

UNIPROT ID:

Q8VCM7

Full Name:

fibrinogen gamma chain

Calculated MW:

49kDa

Purification Method:

Protein A purification

CloneNo.:

241958E11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01455-2, Mouse Fibrinogen gamma chain Recombinant Matched Antibody Pair - PBS only. 84616-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1837. 84616-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL



The mean Fibrinogen gamma chain concentration was determined to be 9,279.6 pg/mL in mouse spleen tissue extract based on a 2.9 mg/mL extract load.