

For Research Use Only

# FADD Recombinant antibody, PBS Only (Capture)

Catalog Number: 84619-1-PBS



## Basic Information

Catalog Number:

84619-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6701

GenBank Accession Number:

BC000334

GeneID (NCBI):

8772

UNIPROT ID:

Q13158

Full Name:

Fas (TNFRSF6)-associated via death domain

Calculated MW:

23 kDa

Purification Method:

Protein A purification

CloneNo.:

241976B8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01426-2, FADD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84619-1-PBS. Detection antibody: 84619-2-PBS. Standard: Ag6701. Range: 0.781-100 ng/mL.